

SEMI AUTOMATIC CAPSULE FILLING



Salient Features:

- cGMP Model
- PLC based machine
- A special design for easy changeover of different capsule size
- Filled weight accuracy within $\pm 3\%$
- Simple operation, operator friendly
- Automatic declutching of loader to minimize operator's attention and ensures optimum production
- Output range upto 30,000 cap/hr




SEMI AUTOMATIC CAPSULE LOADER

Salient Features:

To enhance production line output

- CLM with CFM M/C
- Output: 90 Trays / hour



MANUAL CAPSULE FILLING

Salient Features:

- Output: Upto 7,200 Capsules / Hour.
- Using two operators to load capsule trays and one operator to operate the machine.

ANVAY PHARMA SYSTEMS PVT. LTD.

Manufacturers of Pharmaceutical Process & Allied Machineries



Capsule Section



Liquid Oral Section



Ointment Section



Tablet Granulation Section

Where
Technology
Innovates

ANVAY PHARMA SYSTEMS PVT.LTD.

Unit No. 2 & 3, Novelty Compound, Near Burmashall Petrol Pump, Vasai Highway

Vasai East, Distt. Palghar - 401 208, (Mum.) India

Telephone 91- 250 - 6599551 / 2, Mobile : +91 8888944666

Email : info@anvaypharma.com / sunilchuri@anvaypharma.com

www.anvaypharma.com

www.anvaypharma.com



FLUID BED DRYER

Salient Features:

- Very versatile process equipment
- Fully automatic operation based on PLC control
- Minimum handling and more hygienic
- Top Mounted telescopic cylinder for auto bag shaking
- Clamp type bottom mesh for container
- MOC SS 304 / 316 / 316L
- cGMP with documentation

RAPID MIXER GRANULATOR



MULTIMILL

OCTAGONAL BLENDER



VIBRO SIFTER



BOWL LIFTER

AUTOMATIC OINTMENT CREAM MFG. PLANT

Application

Pharmaceuticals: Creams, Gels, Sterile & Non Sterile Ointments and Biotech products.

Cosmetics: Cold Creams, Hair Dye Cream, Conditioners, Lotions, Gels, Shampoo, Toothpaste.



Salient features:

Main Vessel:

- Choice of various type of anchor design.
- Homogeniser (option) – Top Entry, side Entry,
- Bottom Entry with circulation facility.
- Hydraulic lid lifting for fast & easy cleaning of top dish & vessel.
- VFD for Anchor agitators to vary the RPM.
- VFD for Homogeniser to vary RPM.
- Specially Design drives assembly.

Wax vessel:

- Open top with top Entry agitator Cowl disc type or paddle with Slow Speed Anchor.
- Close top with top/bottom entry agitator Cowl disc type or propeller type stirrer.

Water Vessel:

- Open top with top entry propeller stirrer
- Close top with bottom or top entry propeller.

AUTOMATIC LIQUID/SYRUP MFG. PLANT



Salient Features:

- The complete line is manufactured under CGMP standard
- Consisting of sugar melting vessel, main syrup manufacturing vessel and storage vessel, filter press and pumps.
- Bottom entry stirrers to restrict overall ceiling height.
- Externally placed inline homogenizer ensures the entire material homogenization and gives option for desired multiple passes in recirculation mode with self inbuilt pumping facility within the unit.
- Working platform in S.S. 304. ladder with railing at three sides of the platform will be provided.